



Important Information Regarding Hydrogenics 60

October 23, 2009

Dear Walmart & Sam's Club Independent Doctors of Optometry,

It is our goal to provide you with the highest quality of products and services. With the new advances being made in contact lens technology, we have decided to discontinue the Hydrogenics 60 brand on **October 31, 2009**.

In the coming months we suggest that you refit any of your patients wearing Hydrogenics 60 into one of the following brands that will offer them excellent comfort and health they have come to expect from CooperVision lenses.

Avaira, CooperVision's 2 week silicone hydrogel lens offering your patients:

- Naturally wettable material – no surface treatments or wetting agents
- High Oxygen transmissibility - DK/t 125 helps to maintain healthy corneal physiology
- Low modulus for great comfort

UltraFlex XC, offers CooperVision's proven PC Technology that features a unique biocompatible lens material for all-day and end-of-day comfort. It is also a great option for patients for whom silicone hydrogels are not suitable

UltraFlex II Aspheric, offers your patients excellent clarity of aspheric optics with CooperVision's advanced front surface aspheric design.

As part of our ongoing commitment to the success of your practice, enclosed you will find product specification information on suggested lenses to help you select the best lens for your patients' demands.

To learn more about these products or to request a trial fitting set, please contact your CooperVision Sales Representative or call 800-341-2020. You can also learn more about our products by visiting www.coopervision.com.

Thank you for your business and continued support.

Sincerely,
Sales and Marketing
CooperVision Inc.

MAKE SURE YOUR PATIENTS CONTINUE TO ENJOY THE BEST IN VISION CORRECTION.

To help ensure that your patients currently wearing Hydrogenics® 60 continue to enjoy the high level of vision correction they've come to expect, the chart below offers you other CooperVision lenses that closely match Hydrogenics' parameters and performance.

	MATERIAL	WATER CONTENT	MODALITY	BASE CURVE	DIAMETER	POWER RANGE
HYDROGENICS® 60 UV This product is being discontinued.	Ocufilecon F	60%	Extended wear	3.75 SAG (8.5 mm equivalent)	14.1 mm	-0.25D to -10.00D (0.50D steps over -6.00D)
				3.74 SAG (8.6 mm equivalent)	14.3 mm	+0.25D to +5.00D

We recommend that you fit your patients in one of the lenses below.

Fitting sets and trial lenses are available for your use. Please contact your CooperVision sales representative or authorized distributor for more information.

AVAIRA® with Aquaform® Comfort Science™	Enfilcon A	46%	Daily wear	8.4 mm	14.2 mm	+0.25D to +8.00D (0.50D steps after +6.00D)
				8.5 mm	14.2 mm	-0.25D to -12.00D (0.50D steps after -6.00D)
ULTRAFLEX® XC with PC Technology™	Omafilecon A	60%	Daily wear	8.5 mm	14.2 mm	Plano to -10.00D +0.50D to +6.00D (0.50D steps above -6.00D)
ULTRAFLEX® II ASPHERIC	Ocufilecon D	55%	Conventional extended wear	8.6, 8.9 mm	14.2 mm	-0.25D to -10.00D (0.50D steps above -6.00D)
				8.8 (+) mm	14.2 mm	+0.25D to +6.00D (0.50D steps above +5.00D)

**IF YOU HAVE ANY QUESTIONS, CONTACT YOUR COOPERVISION AREA SALES MANAGER
OR CALL COOPERVISION CUSTOMER SERVICE AT 800-341-2020.**

coopervision.com

800.341.2020



SEE BEYOND THE ORDINARY®